


**Plate with ellipses
Area measuring tool**


Total treatment time = $t_l \times A/ERA$

ULTRASOUND TREATMENT

t_l = local treatment time
 ERA = Effective Radiation Area = transducer size
 A = size of the lesion = skin area



Please note:
 The accuracy of this tool
 may be distorted by
 printing, copying,
 software-modification
 and other processes
 carried out by the users!
 Check out the control
 measures!



© K.Grey.
 VODSKOV, nov.93.

Control
 measures

← 5 cm →

↑ 10 cm ↓